

# CYCLONE BVM

## Most Compact, Easy-to-Use BVM:

The Cyclone® Pocket BVM® is today's most compact, easy-to-use BVM available to provide ventilatory support to patients who are not breathing or who are in respiratory distress.

Its unique collapsible design is quick-to-deploy, easy-to-operate and does not need to be manually inflated prior to use. To help reduce rescuer fatigue, this BVM requires minimum pressure for compression of the bag to ventilate and quickly recovers its shape during use. The bag also has a textured grip to reduce slippage from a gloved hand.

Packaged in an Olive Green rugged, hard plastic, twist-top case, the Cyclone® Pocket BVM® comes complete with oxygen tubing for delivering supplemental oxygen when possible.

The Cyclone® Pocket BVM® is a single use disposable device.

*The compact, fully functional adult BVM that reduces storage space by up to 75% compared to traditional manual resuscitators.*



Rugged olive green hard case protects the device - ideal for the tactical environment ▶

(1) Rugged, twist-top Hard Case with instructions

(1) Fully-functional, Adult Collapsible BVM

Includes Oxygen Inlet, Reservoir Bag and 6 ft 6 in. of Oxygen Tubing

(1) Ergonomic cushioned, see-through, molded silicone Face Mask

**Unique collapsible design packaged 25% smaller than conventional BVMs**




**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION®**

# HOW TO DEPLOY THE IMPROVED POCKET CYCLONE® BVM:



▲ Separate the Face Mask from the folded Resuscitator Bag



◀ Pull both end valves [Front patient valve housing (air/O<sub>2</sub> exit port) and rear air/O<sub>2</sub> inlet valve] outwards in opposite directions while twisting the Patient Valve

◀ Attach the Face Mask to the Patient Valve Housing

## Features:

- Compact adult BVM that reduces storage space by up to 75% of traditional BVM's
- Quick to deploy and does not need to be manually inflated prior to use
- Comes with 6 ft 6 in. of oxygen tubing for delivering supplemental oxygen when available
- Compact and collapsible BVM
- Packaged in a rugged, hard plastic, twist-top case
- Quickly recovers its shape - reducing rescuer fatigue during use.

## Specifications:

- Comes with 6 ft. 6 in. of oxygen tubing
- Packaged - Diameter 5.3 in., Depth 2.9 in.
- Weight: 17.4 oz
- Olive Green Rugged Case

Item#	NSN#
10-0019	6515-01-568-0193



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) 888.689.6277

PRODUCTS WITH A  
**MISSION®**