


INSTRUCTIONS FOR USE


TWO-HANDED APPLICATION




Route the band around the limb, pass the red tip through the slit of the buckle, and position the **CAT** 2-3" above the bleeding site directly to the skin.*

OR


ONE-HANDED APPLICATION




Insert the injured limb through the loop in the band and position the **CAT** 2-3" above the bleeding site directly to the skin.*




Pull band tightly and fasten it back on itself all the way around the limb, but not over the rod clips. Band should be tight enough that tips of three (3) fingers cannot be slid between the band and the limb. If the tips of three (3) fingers slide under band, retighten and re-secure.



Twist the rod until bleeding has stopped.



Secure the rod inside a clip to lock it in place. Check for bleeding and distal pulse. If bleeding is not controlled, or distal pulse is present, consider additional tightening or applying a second **CAT** above and side-by-side to the first. Reassess.




Route the band between the clips and over the rod. Secure rod and band with TIME strap. Record time of application.


*If you cannot be sure or cannot take the additional time to examine where the bleeding is coming from based on the situation, the **CAT** can be effectively applied over clothing as high on the arm or leg as possible. The **CAT** must NOT be applied over solid objects within the clothing. As soon as the situation permits, the injured limb should be evaluated and the **CAT** re-positioned 2"-3" above the injury directly to the skin.




STORING THE CAT IN THE QUICK LAUNCH CONFIGURATION



Pass the red tip through the slit in the buckle. Pull 8" of band through, fold it back and adhere the band to itself.



Flatten the loop formed by the band. Place the buckle in the middle of the flattened band.



Fold the CAT in half placing the buckle at one end. Make sure the TIME strap is open and in ready position; as shown in the illustration.

The CAT is now ready to be placed in your medical kit.

LIMITATION OF LIABILITY:

C-A-T Resources, LLC, its parent company Composite Resources, Inc., its employees, agents, contractors, suppliers, and distributors shall assume no liability for injury or damages arising from the application or use of the Combat Application Tourniquet® (CAT). The user assumes all risk of liability. The CAT should only be used as directed by user's military service component guidelines, EMS authority, or under the supervision of a physician, and in strict compliance with these Instructions for Use.

STANDARD USE GUIDELINES

When applied in accordance with directions, the Combat Application Tourniquet® (CAT) is a safe and effective device for controlling life-threatening extremity bleeding. The CAT has a Maximum Operating Circumference of 35in (88.9cm). The use of any tourniquet for longer than 2 hours may lead to permanent neurological or muscular damage. If you cannot be sure or cannot take the additional time to examine where the bleeding is coming from based on the situation, the CAT can be effectively applied over clothing as high on the arm or leg as possible. The CAT must NOT be applied over solid objects within the clothing. As soon as the situation permits, the injured limb should be evaluated and the CAT re-positioned 2"-3" above the injury directly to the skin. The CAT should be stored in a cool, dry place away from direct sunlight. Re-use of this single use device may introduce risk of infection. Dispose of used CAT in accordance with local regulations for biomedical waste. Please report all serious incidents in relation to the CAT to the manufacturer and appropriate regulatory authority.



NSN 6515-01-521-7976



PN 30-0001

Tactical Black



INST-001 REV02



Exclusive distribution in the U.S. by:
North American Rescue®, LLC
1-888-689-6277
www.NARescue.com
info@NARescue.com



TACTICAL BLACK
NSN 6515-01-521-7976



TRAINER BLUE
NSN 6910-01-560-2972



EMS ORANGE
NSN 6515-01-612-7129

<p>Combat Application Tourniquet® and CAT are trademarks registered in the United States and certain other countries.</p>	
<p>U.S. Patent Nos. 7,842,067; 7,892,253; 8,888,807; 10,016,203; Canada Patent No. 2,569,550; EU Patent Nos. 1,753,344; 2,425,786; Israel Patent No. 179769; Hong Kong Patent No. 1,163,479; Other Patents Pending</p>	<p>Manufactured by C-A-T Resources, LLC 483 Lakeshore Pkwy, Rock Hill, SC 29730 United States of America 803.325.9300</p>