PNEUMERIC™

CAPNOSPOT™

INSTRUCTIONS FOR USE

Device Description

The Pneumeric Capnospot Pneumothorax Decompression Indicator is a capnography device that provides objective, real-time feedback during decompressive thoracostomy. Specifically, the device detects the presence of carbon dioxide in the gas from the decompressed pneumothorax, permitting more accurate placement of the needle or thoracostomy device compared to the current standard of care based upon auditory assessments.

The Pneumeric Capnospot is a therapeutic device intended to guide the placement of pneumothorax decompression needles and thoracostomy devices.

The product consists of a cylindrical polycarbonate tube containing a carbon dioxide detecting paper. The device has a male luer connector on the distal end of the device to connect to decompression needles and thoracostomy devices, a diaphragm to prevent backflow of gases, a color changing indicator and a female luer on the proximal end of the device that can be used to connect accessories.

The product is shipped non-sterile and labeled for single use only. The device is supplied in a single pack configuration labeled with a QR code linking to an electronic Instructions for Use. Each device is sealed in a foil pouch with a desiccant pack.


Figure 1. Pneumeric Capnospot Pneumothorax Decompression Indicator

Indications and Usage

Used for more accurate placement of pneumothorax decompression devices than the current standard of care auditory assessments. The Pneumothorax Decompression Indicator is a colorimetric capnography device indicated to detect the carbon dioxide in the gas passing through the device used to treat pneumothorax with suspected tension physiology.

Contraindications

- Uncooperative patient.
- When pneumothorax decompression is otherwise contraindicated.

Warnings and Precautions:

- Read all instructions carefully. Failure to follow instructions may lead to serious medical consequences.
- Only persons having adequate training in pneumothorax decompression techniques should perform procedure.
- Capnospot has not been tested in patients with open pneumothorax, where treatment should be based on accepted guidelines and standard practices.
- Capnospot should be used in combination with clinical assessments and standard practice guidelines.
- Procedure details will vary according to practitioner's clinical judgement.
- Single Use Only.
- Do not use if package is opened or damaged.
- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.
- Store between 10°F (-12°C) and 140°F (60°C). Storage for extended periods at temperatures below 10°F (-12°C) or above 140°F (60°C) may reduce shelf life.

Method for Use:

Prepare the Capnospot Pneumothorax Decompression Indicator for use:

1. Inspect the pouch and verify that it is unopened and undamaged. **DO NOT** use the device if pouch is found opened or damaged.
2. Gently open the pouch by tearing the notch below the seal and inspect the components for damage. **DO NOT** use the device if damage is evident.
3. The device consists of cylindrical polycarbonate tube with a male luer on the distal end, female luer on the proximal end, and blue indicator paper inside the tube.
4. A minimal amount of green-yellow color in the indicator paper is a normal condition. Do not use the device if majority of the indicator paper is yellow in color.
5. Connect the male luer on the Capnospot to the female luer portion of a needle or thoracostomy device as shown in figure 2 below. **NOTE: CAPNOSPOT MUST BE CONNECTED TO THORACOSTOMY DEVICE BEFORE PROVIDER ATTEMPTS TO DECOMPRESS PNEUMOTHORAX WITH**

SUSPECTED TENSION PHYSIOLOGY. ALTERNATIVELY, IF A NEEDLE ANGIOCATHETER IS BEING USED ONCE THE NEEDLE IS REMOVED, CAPNOSPOT SHOULD BE PROMPTLY ATTACHED TO THE ANGIOCATHETER.


Figure 2. Connection of Decompression Needle/Thoracostomy Device to Capnospot

6. Place the needle or thoracostomy device as clinically indicated for tension pneumothorax.
7. Observe for color change within 10 seconds. The indicator paper should change from blue to yellow if greater than ambient CO₂ levels are detected.
8. If color changes, continue with current clinical care.
9. If no color change occurs, keep the Capnospot connected to the needle or thoracostomy device, evaluate the clinical condition, and re-assess if needle decompression is indicated. Repeat steps 6-8 if required.

Disposal

- The carton is recyclable. Dispose of all packaging materials as appropriate.
- Use solid biohazard waste procedures to discard device.

Symbol Definitions

The following symbols may appear on the device packaging:

Symbol	Definition
	Manufacturer
	Reference Number
	Product Lot Number
	Use By Date (do not use the device after the date shown)
	Do not reuse
	Keep away from sunlight/Keep away from heat
	Caution: Federal (USA) law restricts this device to sale by or on order of a physician.
	Keep Dry
	Consult Instructions for Use
	Temperature Limit
	Do not use if package is damaged

Pneumeric, Inc.

1 (507) 206-1942


823 4th St SW

Patent# US 10,335,524; US 10,842,920 Additional patents pending

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