Emergency Trauma Dressing & Gauze

Instructions for Use:

Intended Use:

- Intended for external hemorrhage control, placed directly on or around a casualty wound.
- Assess casualty wound size, depth, location and contamination. Treat casualty wound with device, wrapping extremity and fasten in place with hook and loop and/or C Clip.

Warning:

- Contents sterile unless packaging open or damaged. If packaging is opened or damaged DO NOT use device.
- Re-use of this device will degrade the efficacy, resulting in adverse casualty reaction, including potential death.
- Continually monitor casualty to ensure bandage is functioning per medical protocols.
- Monitor casualty for tourniquet effect and release bandage pressure as condition warrants.
- Replace the dressing as per medical protocols dictate or <24 hours.
- In the event of a malfunction follow local protocols and report any serious incident to North American Rescue or authorized representative and the competent authority of the Member State.
- Dispose contaminated device in accordance with universal protocols.

Adverse Reaction:

- Failure to adequately apply the bandage may allow for continued bleeding. Overtightening the bandage may cause a tourniquet effect, or tissue ischemia.
- Ensure casualty wound is completely covered and continually monitor casualty per medical protocol.
- Replace bandage as dictated by medical protocols.
- Do not use on casualties with known or suspected allergy to cotton, polyester or spandex.
- Monitor casualties with compromising clinical conditions.
- Evacuate casualty to secondary treatment facility, advising follow on care to monitor casualty for allergic reactions or infections.

Contraindications:

- Within the European Union (EU) the Emergency Trauma Dressing is not indicated for pediatrics or pregnant women due to EU regulatory requirements.
- Device is contraindicated for casualties with allergies to cotton, polyester or spandex.
- Device is not intended for implantation.

MDI Europa GmbH
Langenhagener Str. 71
30855 Hannover-Langenhagen
Germany

Lot number

Consult instructions for use

Manufacturer

Do not resterilize

Do not use if package open or damaged

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