

EVERY DAY CARRY WALLET *EDC*

REV090922

The North American Rescue Every Day Carry (EDC) Wallet was designed to provide users with a versatile, concealable, small carrier for essential point of injury medical products. This wallet fits comfortably in a cargo or pant pocket with a slim cube space and weight. Ideal for low visibility carry of key lifesaving kit contents such as NAR Wound packing gauze, Mini ETD or flat elastic wrap, Hyfin chest seal (compact or regular), flat tape, flat pack or rolled gloves, etc. Large 3" elastic loops on both inside panels allow users the ultimate flexibility in packing options to meet their needs.

Purchasing options include the Wallet as a stand alone item (with no medical equipment) or Kit options of Wound Packing Gauze or Combat Gauze.

Features include Laser Cut Laminate (LCL) 500D nylon construction for low weight and extreme durability. One size fits all, non-reflective medical cross on front, 2" wide Velcro closure to accommodate various packing configurations.

Specifications:

- Durable Laser Cut Laminate nylon construction
- Concealed storage for life-saving medical gear
- Large internal elastic loops for customization
- Slim light weight design

Dimensions:

- Open - L 12.75 in. x H 4.0 in.
- Closed - H 6.25 in. x L 4.0 in. x D 0.35 in.
- Weight: 1.74 oz



Every Day Carry Wallet

ITEM #	COLOR
80-1237	Black
80-1238	Coyote
80-1239	Ranger



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®