

# NAR READY HEAT 4-CELL

REV092922



## 10 hours of continuous dry heat...

The NAR Ready Heat® blanket utilizes proprietary technology and all natural ingredients to deliver continuous and consistent warming. This product is utilizing the same technology carried by every US Army Medic. As an individual unit, The NAR Ready Heat Blanket®'s self-heating, 4-cell liner provides 10 hours of continuous dry heat up to 125°F (52°C). **WARNING: DO NOT PLACE DIRECTLY ON EXPOSED SKIN.**

### Features:

- Blanket begins heating immediately upon opening, no activation required
- Carried by every US Army Medic
- Approved for use on fixed wing and rotary aircraft by the US military
- Self-heating, 4-cell liner
- Provides 10 hours of continuous dry heat up to 125°F (52°C)
- NOTE: DO NOT PLACE DIRECTLY ON EXPOSED SKIN

### Specifications:

- Packaged:  
L 15in. x W 12.75in. x D 1in.
- Blanket:  
L 26in. x W 22in. x D 1in.
- Weight: 1.88 lb oz

### Ready Heat Blanket

| ITEM #  | NSN# |
|---------|------|
| ZZ-1467 | ---- |



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**

WITH

**MISSION®**