

NAR COMPRESSED GAUZE

REV031824



The ultra-compact 100% cotton gauze essential for every medical kit...

North American Rescue's Compressed Gauze is an ultra-compact roll of sterile, high quality 100% cotton gauze that utilizes a unique crinkle weave that has excellent fluid absorption.

Indicated for the treatment of bleeding and the dressing/bandaging of wounds, NAR's Compressed Gauze features a proprietary packaging process that ensures the most consistent, low-cube package on the market.

The package is a rugged and durable vacuum sealed container with easy-to-open tear notches. The ultra-compact storage size expands once opened to cover larger wound areas. NAR's Compressed Gauze is an essential item for any medical kit and can be used to control hemorrhage in conjunction with a compression bandage, used as backing gauze for hemostatic dressings or for wound bandaging.

• Features:

- The compressed, rolled gauze is ideal for wound bandaging and wrapping of injuries
- Unique crinkle weave ensures excellent fluid absorption
- Ultra-compact size that expands to cover larger wound areas
- Used as backing gauze for hemostatic agents, in conjunction with a compression bandage to control hemorrhage

• Specifications:

- 100% Cotton Gauze - 6 ply, W 4.5 in. x L 4.1 yd.
- Packaged: L 3 in. x W 2 in. x D 1 in.
- Weight: 1.2 oz

NAR Compressed Gauze

ITEM #	NSN #
30-0052	6510-01-503-2117



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**