

T3™ LIGHTWEIGHT Assault Litter

REV070822



Features:

- Weighs less than 12 lb
- Maximum maneuverability during TACEVAC in urban terrains
- Overall length of 81 in. enables floor loading on a UH-60 helicopter
- 2 straps for securing patients
- 4 IV attachment points
- Ergonomically designed handles
- Quad-fold design for maximum cube space storage

Our Newest Evolution of Evacuation Platforms..

Weighing less than 12 lb! The T3 Lightweight Assault Litter is designed to support Special Operation Forces where every pound matters.

State of the art materials and technologies have allowed us to reduce the weight without compromising the overall strength and performance of the evacuation platform.

The overall length of 81 in. enables floor loading on a UH-60 helicopter while allowing for maximum maneuverability during casualty evacuation in urban terrain.

Weighing in at less than twelve pounds, the T3™ is the lightest rigid litter in the world.

Specifications:

- Folded: H 21 in. x W 8.75 in. x D 7.25 in.
- Open: L 81 in. x W 22.5 in. x D 6 in.
- Weight: 11 lb 13 oz (5.34 kg)
- Payload: 1200 lb (544.3 kg)
- Ground Clearance: 1.5 in.

T3™ Lightweight

| Item# | NSN# |
|---------|------------------|
| 60-0022 | 6530-01-663-1222 |

Ultra-lightweight, ultra-strong litter weighing less than 12 pounds - specially designed for Special Operations Forces.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®