


STINGRAY® Poleless Litter · Non-Rigid Litters

REV070822



Tested and Approved Poleless Litter...

Constructed with a non-skid, flame-retardant monofilament polypropylene material. Two web straps secure the litter around the body; six heavy-duty handles make it easy to carry.

This design allows rescuers to insert a 2 in. pole, such as pike poles, through the litter cover or to use debris through the outer handles. The Stingray® is ideal for decontamination roller systems and moving wet patients.

Stingray® is tested and approved by both Federal and private laboratories using live warfare agents.

*Resistant to mildew, acids,
alkalis and penetration liquids.*

Specifications:

- Folded: H 20 in. x W 12 in. x D 2 in.
- Open: L 78 in. x W 24 in.
- Weight: 1 lb 13 oz (0.7 kg)
- Payload: 1000 lb (453.6 kg)

StingRay Poleless Litter



Item#	NSN#
60-0004	6530-01-492-4079
60-0004-EU	

CE

*Weighing less than 2 pounds with a
payload of 1000 pounds, the StingRay®
Litter is strong and lightweight.*

Features:

- Constructed of non-skid, flame-retardant monofilament polypropylene fabric
- Resistant to mildew, acids, alkalis and penetration liquids
- Tested by Federal & private labs using live warfare agents
- 6 heavy-duty carrying handles
- 2 casualty securing straps
- Made for decontamination rollers
- Great for long term storage



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®