

CELOX™ -A HEMOSTATIC AGENT

REV071122



CELOX™ -A Hemostatic Granules with Applicator:

CELOX,™ a clinically-proven family of chitosan-based products certified safe-to-use and successful in stopping life threatening hemorrhage on the battlefield, in serious emergencies and in the workplace. From the original granular form to the various CELOX™ gauzes, NAR offers the full spectrum of this trusted hemostatic to meet every medic's need in the pursuit of saving lives from severe trauma.

The unique CELOX™-A applicator is designed to give the user an easy, safe and effective way to stop life-threatening bleeding from difficult to treat penetrating trauma. The instinctive applicator delivery system makes it possible for the CELOX™ Granular hemostatic agent to reach deep into small, penetrating traumatic wounds, getting straight to the source of the bleeding.

Features:

- Stops life-threatening bleeding resulting from shock as found in severe traumatic injuries
- Chitosan-based, works independently of the body's normal clotting mechanisms
- Published success in the field
- Generates no heat

Specifications:

- 6 grams of granules per package
- Pkg: H 8.9 in. x W 4 in. x D 1 in.
- Weight: 1.1 oz.

CELOX™-A Hemostatic Agent

ITEM #	QTY	NSN
30-0067	1	----
30-0140	50	6510-01-623-9934



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®