

NASOPHARYNGEAL AIRWAY

REV042022



Nasopharyngeal Airway with Lubricant

This product shows North American Rescue's attention to details. This NPA was designed to be dependable, efficient, comfortable, and secure.

Dependable. Remains steady during transport. It is not easily dislodged.

Efficient. The product includes a conveniently attached sealed packet of surgical lubricant to facilitate easy insertion.

Comfortable. More compatible with human anatomy than the standard oropharyngeal airway, this product is recommended for use as a nasal airway adjunct in semi-conscious or unconscious patients with an intact gag reflex.

Features:

- Not easily dislodged during transport
- Better tolerated than an oropharyngeal airway by semi-conscious or unconscious patients with an intact gag reflex
- Packaged with surgical lubricant attached
- Multiple sizes available



Specifications:

- Packaged: L 6 in. x W 1.6 in. x D 0.5 in.
- Device Sizes available: 20 French - 36 French
- Weight: 0.5 oz

Item#	Size	Device Length	NSN#
ZZ-1360	20F	100 mm	
ZZ-0424	22F	105 mm	
ZZ-0422	24F	105 mm	
ZZ-0387	26F	110 mm	
ZZ-0034	28F	117 mm	6515-01-529-1187
ZZ-0423	30F	118mm	
ZZ-0388	32F	140 mm	
ZZ-0776	34F	155 mm	
ZZ-0777	36F	170 mm	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®