

# 10 AMP SIMULATOR BATTERY

REV072022

This long-lasting 10-Amp Nickel Metal Hydride (NiMH) battery allows users up to six hours of training on any of the TOMManikin® simulators. The battery will fully charge from a 20% level using the smart charger (not included) in four hours, creating shorter delays in training sessions.

## **Specifications:**

- Long lasting high-output 10-AMP NiMH battery
- Compact and lightweight construction
- Approximate 6 hours of working time on full charge
- Approximate 4 hour full charge from 20% level
- TOMManikin® consumable

## **Dimensions:**

- Battery: L 4 in. x H 6 in. x W 6 in.
- Weight: 4.0 lbs

## **10 Amp Simulator Battery**

ITEM #	NSN #
93-0141	--



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**