

TOMMANIKIN® GSW with Amputation

REV092322

Enhances TCCC Instruction

TOMManikin® GSW (Gun Shot Wound) with Amputation features through-through GSW to leg above the knee. Entry wound presents with small entrance hole with circular powder burn. Exit wound has exposed fat and wound site ecchymosis. Variable bleeding rate. Responds to both peripheral and junctional tourniquet application and femoral pressure point application.

Leg amputation presents with peripheral wound site ecchymosis, exposed sub-q fat layer and muscle striation. Entry wound is packable and responds to both peripheral and junctional tourniquet application and femoral pressure point application

All TOMManikins include: NARS Tablet (software includes: evaluation checklists, scenario presets and scenario builder software), extended care operation kit, 1-year warranty, TOMManikin 5 Amp and 10 Amp batteries with smart chargers, onboard fluid fill unit assembly tank, 8 oz simulated blood concentrate, 5 surgical airway neck skins (35 op sites total), 2 3-wheeled custom deployment bags for safe storage and transport, Velcro uniform, NAR IFAK.

Features:

- Industry-leading durability - withstands drops/falls/drag
- Easy set-up, rapid reset and low maintenance
- Sturdy steel frame and joints
- Weighs between 150-170 lbs. (depending on configuration)
- Water Resistance Tactical Operation Manikin & remote-control operating system
- Rugged Softech®, non-leaching skin and extremities
- Massive bleeding from multiple locations
- Realistic sucking chest wound (standard on GSW platform)
- Quick replacement of extremities for flexibility with multiple wound configurations
- Interchangeable extremities between TOMM platforms (GSW, Blast, Basic, TAMIKIN)

- Manikin Voice with the ability to upload multiple custom voice tracks in any language
- Voice tracks can be remotely changed on the fly. Onboard SD card stores up to 16 GB data
- Eyes are interchangeable to reflect changing pupil conditions
- Onboard battery is a 12v / 5 amp & 10-amp standard nickel metal hydride (NiMH) batteries
- Bends and rotates for realistic patient movement at the arms, legs and waist
- Realistic bilateral rise and fall of chest
- Pelvic crepitus
- Remote control allows instructor to control talking, breathing and bleeding and response to interventions
- Features proprietary communication link that requires NO Wi-Fi, Bluetooth, or LAN connections
- Simultaneous long-range communications with numerous NARS devices

Dimensions:

- L 40 in. x H 36 in. x W 48 in.
- Weight: 240 lbs

TOMMANIKIN® - GSW with Amputation

ITEM #	STYLE	NSN #
93-0369	GSW w/ AMP	6910-01-705-3130



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277



PRODUCTS
WITH A
MISSION®

TOMMANIKIN® GSW with Amputation

REV092322

TCCC Interventions Supported:

- Massive hemorrhage control (tourniquet application and wound packing)
- Chest seal application
- Pleural needle decompression (multiple sites)
- Airway management NP/OP airways, supraglottic devices, endotracheal intubation, surgical airway, and chest tube (Upon request)
- Tibial intraosseous site

Kit Contents:

- 1 x Remote instructor tablet
- 1 x Blood reservoir pressure fill tank
- 1 x Velcro digital pattern uniform
- 2 x 3-wheeled deployment-style bags for safe storage and transport
- 1 x NAR Training IFAK
- 1 x Operator's Manual

Injuries Included:

- Features anterior GSW entry & mid-auxiliary exit. Actively bleeds from both locations. Anterior entry wound mixes with air from onboard compressor to create bubbling bleeding effect. Supports chest seal application and pleural needle decompression
- TOMManikin GSW chest skin
- TOMManikin GSW pelvic skin
- Features GSW anterior entry/posterior exit with soft tissue trauma to penis and avulsed scrotum with prolapsed testicle
- Leg amputation present with peripheral wound site ecchymosis, exposed sub-q fat layer and muscle striation. Entry wound is packable and responds to both peripheral and junctional tourniquet application and femoral pressure point application. This wound is interchangeable with other TOMManikins
- The TOMManikin Uninjured Arms feature radial pulse, antecubital IV access and humeral IO site
- Unstable pelvis with palpable crepitus



TOMMANIKIN® - GSW with Amputation

ITEM #	STYLE	NSN #
93-0369	GSW w/ AMP	6910-01-705-3130



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®