

# BLACK TALON® GLOVES

REV120222



*Sold on bulk box  
or in pre-packaged  
pairs of 25*

## High performance TACTICAL BLACK textured gloves for improved grip in wet and dry situations

Black Talon® Ultimate Nitrile gloves are made of 100% nitrile material which is free of the allergens found in standard latex gloves. They are engineered to give our users high performance by offering both maximum protection and unsurpassed tactile sensitivity.

This careful attention to detail includes textured fingertips to help with fine motor skills, beaded cuffs for quick and easy donning and a thickened palm for durability and reliability. These ambidextrous gloves are packaged 50 pairs per box (Extra large, 45 pairs per box); 25 pairs per bag. Designed for single-use only.

### Features:

- 100% nitrile gloves, textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Extended, beaded cuff for quick and easy donning
- Available in Tactical Black for reduced exposure when stealth is required
- Form fitting, Non-sterile, Available in sizes S, M, L and XL



### Glove Specs:

- Length: 11.8 in.
- Cuff Thickness (mils): 4.7
- Palm Thickness (mils): 6.0
- Finger Thickness (mils): 8.0

### Bulk Dimensions:

- (Box) L 11 in. x W 5 in. x D 3.75 in.
- Weight: 2 lb 4 oz

### Pack Dimensions:

- (Bag) L 4.75 in. x W 6 in. x D 4 in.
- Weight: 16.1 oz

Size	Unit	Qty.	Item#	NSN#
S	Box	50 pr.	70-0001	6515-01-515-0195
S	Pack	25 pr.	ZZ-0218	
M	Box	50 pr.	70-0002	6515-01-515-0200
M	Pack	25 pr.	ZZ-0217	6515-01-672-7977*
L	Box	50 pr.	70-0003	6515-01-515-0197
L	Pack	25 pr.	ZZ-0216	6515-01-632-0059
XL	Box	45 pr.	70-0004	6515-01-515-0190
XL	Pack	25 pr.	ZZ-0215	6515-01-672-7985*

\*NATO NSN



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**