

TTMB Tiered Trauma Medical Bag

REV112717

5-Piece Modular Ensemble

The five-piece modular TTMB (Tiered Trauma Medical Bag) targets the military medic's requisite that offers countless load carriage platform options for mission specific applications.

This ensemble offers load carriage platform options for mission specific applications, including a low profile backpack style aid bag, a re-supply bag, a medical chest pouch, and a leg rig conversion panel with a MOLLE directional conversion panel. The ample re-supply bag includes modular medical supply pouches that can be used to restock both the aid bag and chest pouch. The leg rig conversion and MOLLE directional conversion panels allow the chest pouch to be mounted and carried using a drop-down leg attachment style platform.



**Nylon
Colors:**



Kit Contents:

- 1 x ReSupply Bag
- 1 x Aid Bag
- 1 x Chest Pouch
- 1 x Leg Rig Conversion Panel
- 1 x MOLLE Directional Conversion Panel

NOTE: Medical equipment is not included.

ORDERING INFORMATION

BAG	Color	ITEM #
TTMB Complete System	Multi-Cam	80-0274
ReSupply Bag (A)	Multi-Cam	80-0271
NAR-4 Aid Bag (B)	Multi-Cam	80-0241
NAR-4 Chest Pouch (C)	Multi-Cam	80-0248
Leg Rig Panel (D)	Multi-Cam	80-0272
MOLLE Conv. Panel (E)	Multi-Cam	80-0273

NOTE: Bags only all contents sold separately

Features:

- 5 piece modular system with inter-changable inner pouches compatible with any Hook-&-Loop attachment system
- Removable labels aid in the organization of mission critical medical equipment
- Constructed of 500D CORDURA® Berry Amendment Nylon, solution dyed for near IR signature reduction
- GhillieTEX hardware, sand & water resistant self-repairing zippers & silent 550 cord pulls
- Includes low-profile backpackstyle aid bag, ample re-supply bag, medical chest pouch, leg rig conversion and MOLLE directional conversion panels
- Bags only. Medical equipment is not included

Dimensions:

- ReSupply Bag: H 20 in. x W 15 in. x D 5 in.
- Aid Bag: H 19 in. x W 15 in. x D 7 in.
- Chest Pouch: H 7.75 in. x W 9.75 in. x D 4.5 in
- Leg Rig Conversion Panel: L 14 in. x W 6 in.
- MOLLE Directional Conversion Panel: L 7.5 in. x W 6 in.
- Total Weight: 9 lb 2 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

TTMB Tiered Trauma Medical Bag



A. ReSupply Bag

Item# 80-0271
(MTC only)

C. NAR-4 Chest Pouch

Item# 80-0248 (MTC)
Other colors available

B. NAR-4 Aid Bag
Item# 80-0241 (MTC)

Other colors available
see NARescue.com

SPECIFICATIONS

- (A) ReSupply Bag: H 20 in. x W 15 in. x D 5 in.
- (B) NAR-4 Aid Bag: H 20 in. x W 15 in. x D 5 in.
- (C) NAR-4 Chest Pouch: H 7.75 in. x W 9.75 in. x D 4.5 in.
- (D) Leg Rig Conversion Panel: L 14 in. x W 6 in.
- (E) MOLLE Conversion Panel: L 7.5 in. x W 6 in.



D. Leg Rig Conversion Panel
Item# 80-0272 (MTC only)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277