

# K-9 HANDLER IFAK

REV112624



## ***Designed for Immediate Lifesaving Treatment to a Wounded Handler, Operator or K-9 at Point of Injury...***

The NAR K-9 Handler Kit is designed to enable canine handlers to provide lifesaving treatment to injured canines at the point of injury or when veterinary care is not immediately available. This kit and its contents are streamlined, durable and allow the handler to carry an individual first aid kit that can be utilized for rendering care to the handler, their canine or other injured persons.

The NAR K-9 Handler Kit provides the necessary tools to address the major traumatic injuries seen in working dogs (including the option of choosing the hemostatic dressing Combat Gauze®). It also includes treatment for other canine emergencies such as bloat or heat injury. Small size allows for easy attachment using MOLLE/PALS style connectors to the assault vest/body armor or can attach to a Quick Detach Panel for fast, easy removal of entire kit (available for separate purchase).

### **Features:**

- Rugged IFAK for point of injury care to canines and/or humans
- Address canine emergencies from major bleeding to bloat or heat injury
- MOLLE/PALS to assault vest/body armor for easy attachment
- Includes the fastest, safest and most effective prehospital tourniquet in the world, the C-A-T® – for treating handlers or injured humans

### **Specifications:**

- Kit Closed:  
H 7 in. x W 4.5 in. x D 2.75 in.
- Weight: 1 lb 5 oz

### **K-9 Handler IFAK**

COLOR	OPTION	ITEM#	NSN#
BLK	Basic	80-0439	--
COY	Basic	80-0440	--
MTC	Basic	80-0441	--
ODG	Basic	80-0442	--
RGR	Basic	80-1504	--
RTC	Basic	80-1505	--
RED	Basic	80-1503	--
BLK	Advanced	80-0443	--
COY	Advanced	80-0444	--
MTC	Advanced	80-0445	--
ODG	Advanced	80-0446	--
RGR	Advanced	80-1507	--
RTC	Advanced	80-1508	--
RED	Advanced	80-1506	--



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# K-9 HANDLER IFAK

REV112624



**+plus**



*Address Canine Emergencies from Major Bleeding to Bloat to Heat Injuries.*

## Kit Contents:

- 1x C-A-T® Tourniquet- BLK
  - 1x 4 in. Emergency Trauma Dressing (ETD)
  - 1x HyFin® Vent Chest Seal, Twin Pack
  - 2x Petrolatum Gauze
  - 8x Alcohol Prep Pad
  - 2x Pair, Bear Claw Nitrile Gloves, Lg
  - 1x Surgical Tape, 1 in.
  - 1x Digital Vet Thermometer
- Basic kits also include:**
- 2x NAR Compressed Gauze
- Advanced kits also include:**
- 1x NAR Compressed Gauze
  - 1x Z-Fold Combat Gauze®
  - 1x ARS®, Needle Decompression Kit, 14 ga.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A **MISSION®**