

X-SHEARS® SOFT HOLSTER

REV021221



Shears sold separately



Carry your Shears in comfort...

The X-Shears® Soft Holster is designed for comfort and easy access. Durable fabric over soft foam makes this holster comfortable no matter how you wear it. The powder coated heavy steel belt clip will grip even thin fabrics and may be worn inside or outside of pockets, providing ready access to the shears while protecting your clothing.

The XShear Soft Holster is designed with XShear Trauma Shears in mind but compatible with many types of shears.

Cleaning - Hand wash with cleaning solution or disinfectant. Rinse and allow to air dry.

Outside material is Spandex with a foam core and nylon inside material.

Features:

- Designed for comfort
- Convenient solution for easy access to shears
- Allows Left or Right hand carriage
- Powder coated heavy steel belt clip

Specifications:

- Dimensions:
L 7 in. x W 3.5 in. x H 0.5 in.
- Weight: 2 oz

X-Shears Soft Holster

ITEM #	NSN #
70-0314	---



Compatible with many Trauma Shears Sold Separately



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®