

DISPOSABLE SURGICAL MASK

LEVEL 3 (Box of 50)

REV101220



A Must Have...

These disposable surgical masks are a level 3, high barrier mask. With Bacterial Filtration Efficiency (*BFE*) and Particle Filtration Efficiency (*PFE*) at 98% with 0.1 micron, you are assured a high level of protection against potential sources of infection. The pleated design will prevent the pooling of fluids and cross contamination, while dye free non-woven inner layer provides additional protection. Soft and strong, ample-length earloops help provide a comfortable fit.

Features:

- Pleated with elastic earloops for a comfortable fit
- Bacterial Filtration Efficiency (*BFE*) and Particle Filtration Efficiency (*PFE*) at 98% with 1.0 micron
- ASTM Level 3 Protection
- Metal nose clip on top of mask for snug fit across nose
- Convenient alternative to cotton face coverings
- Latex free, single use
- One size: Adult
- Sold as a box of 50

Specifications:

- Dimensions (box of 50): 7.5 in. x 4 in. x 5.25 in.
- Weight (box of 50): 0.61 lb

Disposable Surgical Mask - Level 3

Item #	NSN #
ZZ-1352	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®