

DISPOSABLE SURGICAL MASKS - LEVEL 2/3 *Box of 50*

REV050322

A convenient alternative to cotton face coverings, NAR's disposable Surgical Masks sold in a box of 50 feature ASTM Level 2/3 Protection in a pleated, 3-ply design with elastic earloops. The 3 layers of soft non-woven fabric make these disposable face masks breathable and comfortable to use every day. The metal nose clip on top of mask easily adjusts for a snug fit across nose.

Specifications:

- Pleated, 3-ply design with elastic earloops
- 3 layers of soft non-woven fabric making it breathable and comfortable to use every day
- Metal nose clip on top of mask for snug fit across nose
- Convenient alternative to cotton face coverings
- ASTM Level 2/3 Protection
- Non Sterile, Latex Free, Single Use
- One size: Adult
- Sold in box of 50

Dimensions:

- Mask: L 6.89 in. x W 3.72 in.
- Box of 50: L 7.5 in. x W 4.25 in. x D 4 in.
- Weight: 0.5 lb



Disposable Surgical Masks - Level 2/3

ITEM #	QUANTITY	NSN #
ZZ-1325	Box of 50	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®