

CLASS A & B Kit Hard Cases- ANSI

A portable solution to store & carry First Aid equipment in the workplace or at a work site

This rugged hard case is easily transported to the patient's side and protects contents from the elements, dust and insects. It closes securely using latches, utilizes a gasket to securely seal it, and includes a pressure release valve to ensure that it will easily open if moved between high and low pressure environments. The case offers two locking points that can be used to secure it, or employ tamper seals to ensure contents can be monitored for use. Inner nylon pouches with zipper closures are held in position using hook and loop. Each pouch has a hook and loop patch on the lid for easy identification of the contents. The larger Class B case has a lid that offers clear plastic pouches for easily viewable storage. These cases offer users the ability to customize the position and contents to suit their specific need.

Features:

- Modular Type IV ANSI Container featuring a rugged storage organizer in the lid and removable nylon pouches with easy-to-read treatment labels
- Rugged hard case that protects your medical equipment from the elements
- Designed to be easily carried to the patient's side
- Features pressure equalization valve for air transport
- Can be locked or tamper sealed

Class A & B Kit Hard Cases - ANSI

SIZE	ITEM #	NSN #
Class A	80-1113	
Class B	80-1114	

Dimensions:

Class A

- L 13 in. x H 9.5 in. x D 5 in.
- Weight: 6 lb

Class B

- L 13.5 in. x H 11 in. x D 6 in.
- Weight: 8.87 lb



Class A Case

Class B Case

ANSI
Type IV ANSI Containers



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH

MISSION®