

# AER INFLATABLE CERVICAL COLLAR

REV012723



## ***Latest Technology in Cervical Collars***

The AER Cervical Collar represents the latest technology in cervical spinal immobilization for austere environments, offering lighter weight, smaller packability and greater performance when compared to leading industry competitors.

## ***Specifications:***

- AER collar rigidity can be adjusted by inflating/releasing air
- The AER collar is designed to provide adequate access to the patient for the medical practitioner in all airway emergencies. Quick deflation, adjustment of rigidity or access to the anterior neck is all possible with the AER collar
- Packable size allows for easy carry and storage (2.4 oz)
- Chin Strap provides intuitive placing and comfort of collar
- Easy and secure latching allows for varying sizes
- Air battens create structure while maintaining critical blood flow to the brain
- Durable 75d Ripstop Nylon is strong for most situations and environments

## ***Dimensions:***

- Packaged - L 7.5 in. x H 3 in.
- Weight: 2.4 oz

## ***AER INFLATABLE CERVICAL COLLAR***

ITEM #	QUANTITY	NSN #
50-1010	1 Each	--



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**