

SHIRT SHIELDS for Rigid Gen7 C-A-T® Case

REV061818

Prevents wearing of uniform shirts by covering the windlass clip end...

The RIGID CASE SHIRT SHIELD is an add on accessory that easily mounts to your RIGID TQ CASE. The SHIRT SHIELD attaches easily with two screws and is designed to help prevent wearing of uniform shirts against the windlass clip on your C-A-T® Tourniquet.

Features:

- Prevents wearing of uniform shirts by covering the windlass clip ends
- Easily attaches to any existing RIGID C-A-T® TOURNIQUET CASE with two screws
- SHIRT SHIELD is made from injection molded nylon
- Allows for easy and quick access to your C-A-T® Tourniquet while providing protection for your uniform shirt

Dimensions:

- H 1.75 in. x W 2.5 in. x D 1.5 in.
- Weight: 1.44 oz

Shirt Shield for Rigid Gen7 C-A-T® TQ Case

Color	Item #	NSN#
BLK	80-0941	
COY	80-0942	
RGR	80-0943	



Rigid Gen7 C-A-T®
TQ Case (sold separately)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**