

# RIGID C-A-T® TQ CASE COVER

REV113020

## Helps Protect Your C-A-T® Tourniquet from the Elements...

The C-A-T® RIGID TQ Case Cover is an optional accessory for your C-A-T® RIGID TQ Case. The Case Cover completely covers your tourniquet to protect it from the elements and provides a clean, uniform look.

The cover integrates with the existing hardware on your C-A-T® RIGID's mounting hardware and includes a snap that utilizes an existing hole on the front of the case. When deployed, the cover "breaks apart" so there is no snagging on your tourniquet..

*Please Note: Accommodates both Gen 6 and Gen 7 C-A-T® Tourniquets (not included)*

*Does not include the C-A-T® RIGID TQ Case pictured.*

*Will not work with the Slick Front RIGID TQ Case or the legacy Kydex TQ Cases*


### Features:

- Lightweight Hypalon One-Piece Construction adds little weight to your C-A-T® RIGID TQ Case
- Includes Loop square for a medical identifier patch (not included)
- Contains all necessary hardware to integrate into the C-A-T® RIGID TQ Case with no modifications
- Completely covers the tourniquet
- Break-Away design allows smooth removal of your tourniquet with no snagging
- Accommodates both Gen 6 and Gen 7 C-A-T® Tourniquets (sold separately)

*NOTE: Tourniquets and C-A-T® RIGID TQ Case sold separately*

### Dimensions:

- L 2.5 in. x W 2.125 in. x D 1.5 in.
- Weight: 1 oz



### Rigid C-A-T® TQ Cover

COLOR	ITEM #
Black	ZZ-0661
Coyote	ZZ-0662
Ranger Green	ZZ-0663
MultiCam	ZZ-1366



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**