

Rigid Gen 7 C-A-T TQ Case MOUNTING SCREWS for SAFARILAND Drop Leg Shroud (Set)

REV021221



Keep your Gen 7 C-A-T® Tourniquet accessible!

This item provides the mounting hardware needed to attach an Eleven 10 RIGID TQ Case® to a Safariland® Drop Leg Shroud in the accessory position.

Package includes 2 barrel nuts, bushings, screws and washers. (TQ Holder and Safariland Drop Leg Shroud not included.)

Features:

- Mount your Rigid Tourniquet Holder to the Safariland Drop Leg Shroud
- Includes 2 barrel nuts, bushings, screws and washers

Specifications:

- Packaged: L 5.7 in. x W 4.0 in. x D 0.5 in.
- Weight: 4 oz

Mounting Screws

ITEM #	NSN #
ZZ-0562	---



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®