

CAT[®] ELASTIC TQ HOLDER

REV080718

Rapid, easy access from your vest, gear or belt...

The C-A-T Elastic TQ Holder from North American Rescue was specifically designed to allow personnel to place their C-A-T[®] Tourniquet on their vest, gear or belt for rapid and easy access. It is made of 500D IR Signature Reduced Nylon with MOLLE/PALS-style connectors with end tuck instead of button snap. The Elastic Holder protects your tourniquet from the elements and has an easy-open hood design to completely clear the TQ for removal up down or outward.

NOTE: Accommodates both Gen 6 and Gen 7 C-A-T[®] Tourniquets (not included)

Features:

- Tailored specifically for your C-A-T[®], to protect it from exposure to the elements
- PALS/MOLLE-style strap with tuck finish eliminating the need for a button snap
- Versatile attachment to the duty vest, gear or belt for rapid & easy access
- Easy-open quick release hood
- Extended elastic body for TQ release outwards under pressure
- 500D IR Signature Reduced Nylon
- Accommodates both Gen 6 and Gen 7 C-A-T[®] Tourniquets (not included)

Dimensions:

- L 7 in. x W 1.75 in. x D 1.5 in.
- Weight: 2 oz

C-A-T[®] Elastic TQ Holder

Color	Item #	NSN#
BLK	30-0184	
COY	30-0185	
MTC	30-0186	



NOTE:
Tourniquets
sold separately



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION[®]