

ANKLE TRAUMA KITS (ATK)

REV020518



Keep essential trauma medical equipment readily accessible without adding bulk on the belt...

The Ankle Trauma Kit (ATK) from North American Rescue offers Law enforcement officers or other individuals the capability of carrying their C-A-T® Tourniquet, as well as an additional ETD™ and gauze/hemostatic all in an ultralight ankle holster.

Designed with a lightweight, water resistant ripstop nylon frame and pouches that keep the items dry and secure, the ATK also features a soft, breathable fabric that helps with moisture wicking on the interior for extended user comfort. The ankle holster and trauma kits are “one size fits most”, implementing high strength hook and loop closures for maximum adjustability and versatility. The contents of the ankle holster can be accessed one-handed, with either hand quickly and easily for rapid deployment.

The Trauma Ankle Holster as a stand alone product is also available for separate purchase.



Rapid, single-handed, access to components with either hand...



Equipment Option	Item#	NSN#
Ankle Trauma Kit with Combat Gauze LE	80-0870	
Ankle Trauma Kit with S-rolled Gauze	80-0869	
Ankle Trauma Holster only	30-0177T	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

ANKLE TRAUMA KITS (ATK)



Ankle Trauma Kit with Combat Gauze LE Contents:

- 1 x C-A-T® Combat Application Tourniquet, Black
- 1 x Emergency Trauma Dressing – 4 in. Flat
- 1 x Combat Gauze® LE Hemostatic Dressing
- 1 x Ankle Trauma Holster, Black



Ankle Trauma Kit with S-rolled Gauze Contents:

- 1 x C-A-T® Combat Application Tourniquet, Black
- 1 x Emergency Trauma Dressing – 4 in. Flat
- 1 x S-rolled Gauze
- 1 x Ankle Trauma Holster, Black

Features:

- Very lightweight, low profile, versatile carrier that allows for carry of essential trauma equipment
- Designed to accommodate a C-A-T® Tourniquet, Emergency Trauma Dressing, and Hemostatic/Gauze
- Advanced breathable materials that provide user comfort and long-term wearability
- Protective, water resistant pouches for each carried component with strong hook and loop closures
- Quick, single-handed, ambidextrous access to components
- One size fits most with convenient hook and loop adjustability

Specifications:

- KIT with Combat Gauze LE:
L 17.65 in. x H 6.5 in. x D 1.5 in. | Weight: 8.25 oz
- KIT with S-rolled Gauze:
L 17.65 in. x H 6.5 in. x D 1.5 in.. | Weight: 8.75 oz
- Ankle Holster only:
L 17.75 in. x H 4 in. x D 0.25 in. | Weight: 2 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®