

SPORTSMAN M-FAK *Mini First Aid Kit - REALTREE®*

REV092721

NAR's most compact, versatile Mini First Aid Kit

The Sportsman version of North American Rescue's M-FAK Mini First Aid Kit, was designed to be the most compact, multi-use IFAK for delivering immediate critical care for penetrating or other traumatic injuries encountered in the field. At half the size of our popular T.O.R.K. platform, this M-FAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by First Responders operating in the line of duty.

This M-FAK's super compact, 500 Denier nylon bag comes in REALTREE® Camo or opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T Tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips to attach additional gear such as a TQ/Holder (sold separately) if needed.



Features:

- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating or other traumatic injuries in the field
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access
- Multiple elastic loops for secure gear organization
- Vertical mount can be set to open left-to-right or right-to-left based on shooter preference
- Horizontal mount on a belt allows opening directly to your C-A-T®
- Half the size of our popular T.O.R.K. platform

Available Color:



RTC

Specifications:

- L 6 in. x W 3 in. x D 3 in.
- Weight: 13 oz

Sportsman M-FAK Mini First Aid Kit

| Level | Item# | NSN# |
|--------------------------------|---------|------|
| Basic | 80-0913 | |
| with Bleeding Control Dressing | 80-0914 | |



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

SPORTSMAN M-FAK *Mini First Aid Kit - REALTREE®*



REV092721

**The First Aid Kit
offered in REALTREE® Camo...
Ideal for any sportsman & hunter!**

BASIC -or-
BASIC with Bleeding Control Dressing:



-OR-



The Sportsman REALTREE® Camo M-FAK also includes a PVC MED ID Cross patch



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A  **MISSION®**