

BLACK RESPONDER GLOVES

REV100623



Third party tested for effective protection from Fentanyl as well as bacteria, fungi and viruses including COVID-19 Coronavirus!

Crafted with revolutionary Eco Best Technology® (EBT), NAR's Black Responder Gloves are biodegradable nitrile disposable gloves that deliver unmatched comfort, dexterity and performance across a wide range of applications. These powder free, ambidextrous gloves come 100 gloves per box delivering ultimate PPE performance - third-party tested compliant with the EN ISO 374-5:2016 Standard for Virus Protection against bacteria, fungi and viruses.

Features:

- 100 gloves boxed for rapid deployment (**Note:** XXL has 90).
- High performance personal protection compliant with the EN ISO 374-5:2016 Standard for Virus Protection
- Third-party tested for effective personal protection from Fentanyl and Viruses including COVID-19 Coronavirus
- Latex- and Powder- Free with decreased risk of allergies
- Biodegradable and disposable with second-skin fit & feel
- Black color
- Available in 6 sizes: XS, S, M, L, XL and XXL

Glove Specs:

- Length: 9.5 in.
- Cuff Thickness (mils): 4.6
- Palm Thickness (mils) 4.0
- Finger Thickness (mils): 8.0
- Cuff: Rolled

Box Dimensions:

- L 11 in. x W 4.5 in. x D 3.75 in.
- Weight: 2 lb 5 oz

Black Responder Gloves

Size	Qty.	Item#
XS	50 pr	70-0319
S	50 pr	70-0320
M	50 pr	70-0321
L	50 pr	70-0322
XL	50 pr	70-0323
XXL	45 pr	70-0324



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**