

BLACK RESPONDER GLOVES PACK OF 25

REV100623

COLOR:
Tactical
Black

Third-party tested compliant with the EN ISO 374-5:2016 Standard for **Virus Protection against bacteria, fungi and viruses** including, but not limited to:

- ▶ **COVID-19** (Coronavirus)
- ▶ **HIV** virus
- ▶ **Ebola** virus
- ▶ **Avian influenza** virus (Bird Flu)
- ▶ **Hantavirus** [causes Hantavirus Pulmonary Syndrome (HPS)]
 - ▶ **Marburg** virus
 - ▶ **Rabies** virus
 - ▶ **Smallpox** virus
- ▶ **Influenza** viruses (A, B and C)
 - ▶ **Rotavirus**
 - ▶ **HPV** virus

Third party tested for effective protection from Fentanyl as well as bacteria, fungi and viruses including COVID-19 Coronavirus!

The Medical Protective Face Shield is designed as a cost effective general splash guard to add an extra layer of personal protection from debris and chemical hazards. This professional face shield is simply designed to work in conjunction with a face mask to provide another level of PPE. The OSHA compliant clear shield is constructed of 14 mil BPA free plastic with a fully adjustable, one size fits all headband that can accommodate extra room between the face and the shield for goggles. Made in the USA.

Features:

- 25 individually wrapped pairs for convenient rapid deployment
- High performance personal protection compliant with the EN ISO 374-5:2016 Standard for Virus Protection
- Third-party tested for effective personal protection from Fentanyl and Viruses including COVID-19 Coronavirus
- Latex- and Powder- Free with decreased risk of allergies
- Biodegradable and disposable with second-skin fit & feel
- Available in 6 sizes: XS, S, M, L, XL, & XXL

Specifications:

- Length: 9.5 in.
- Cuff Thickness (mils): 4.6
- Palm Thickness (mils) 4.0
- Finger Thickness (mils) 8.0
- Cuff: Rolled
- Packaged (25 pairs): L 4.75 in. x W 6 in. x D 4 in.
- Individual pair: L 2.5 in. x W 1 in. x D .8 in.
- Weight: 16.1 oz

Black Responder Gloves - Pack of 25

| SIZE | ITEM # |
|------|---------|
| XS | ZZ-1285 |
| S | ZZ-1286 |
| M | ZZ-1287 |
| L | ZZ-1288 |
| XL | ZZ-1289 |
| XXL | ZZ-1290 |



Personal Protection from Fentanyl and COVID-19



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**