


# What Size Do I Wear?

## Sizing Guide:

The size you should choose is indicated by the width of your palm.



**NORTH AMERICAN RESCUE**  
www.NARescue.com • 888.689.6277



**NORTH AMERICAN RESCUE®**


[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# Glove Sizing Chart

## **DIRECTIONS:**

1. Place thumb and forefinger in pattern shown.
2. Note color at tip of finger.
3. Match color to following size category:

■ Yellow	XXS	(05)
■ Green	XS	(06)
□ White	S	(07)
■ Blue	M	(08)
■ Red	L	(09)
■ Black	XL	(10)





**NORTH AMERICAN RESCUE**  
www.NARescue.com • 888.689.6277

## GLOVE SIZING CHART



SIZE	
SMALL	7-8 IN
MEDIUM	8-9 IN
LARGE	9-10 IN
X-LARGE	10+IN



Measure your dominant hand around the fullest part of your hand.  
Do not measure around the thumb.