

RESPONDER NITRILE GLOVES *Blue*

REV083122

Third party tested for effective protection from Fentanyl as well as bacteria, fungi, and viruses, including COVID-19 Coronavirus!

NAR's Blue Responder Gloves are biodegradable nitrile disposable gloves that deliver unmatched comfort, dexterity, and performance across a wide range of applications. These powder-free, ambidextrous gloves come 100 gloves per box delivering ultimate PPE performance—compliant with EN ISO 374-5:2016 Standard for Virus Protection against bacteria, fungi, and viruses, including:

- COVID-19 (Coronavirus)
- HIV virus
- Ebola virus
- Avian influenza virus (Bird Flu)
- Hantavirus
- Marburg virus
- Rabies virus
- Smallpox virus
- Influenza viruses A, B, & C
- Rotavirus
- HPV virus

Specifications:

- 100 gloves boxed for rapid deployment
- High performance personal protection compliant with the EN ISO 374-5:2016 Standard for Virus Protection
- Third party tested for effective personal protection from Fentanyl and Viruses including COVID-19
- Coronavirus
- Latex and Powder Free with decreased risk of allergies
- Biodegradable and disposable with second-skin fit and feel
- Blue color
- Available in 6 sizes: XS, S, M, L, XL, and XXL

Glove Specifications:

- Length: 9.5 in.
- Cuff Thickness (mils): 4.6
- Palm Thickness (mils): 4.0
- Finger Thickness (mils): 8.0
- Cuff: Rolled

Dimensions:

- Packaged: L 4.5 in. x W 11 in. x D 3.75 in.
- Weight: 2.5 lbs



Responder Nitrile Gloves - Blue



SIZE	UNIT	QTY	ITEM#	NSN#
S	Box	50 pr.	70-0315	--
M	Box	50 pr.	70-0316	--
L	Box	50 pr.	70-0317	--
XL	Box	50 pr.	70-0318	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**