

# RAPID RESPONSE KIT *Rescue Task Force edition*

REV100421



BLK

RED

Rapid Response Kits can be modified upon special request. Contact a North American Rescue sales rep for more details.

## **Rescue Task Force Edition is equipped for Initial Triage of a Multiple Casualty Incident...**

Exterior of the rugged nylon shoulder bag features adjustable, padded strap with quick release buckle. Multiple storage capabilities including loop-tape/MOLLE attachment system, 2 extra pockets & elastic storage loops with color-coded labels to accommodate up to 8 chem lights (not included).

Interior compartment features 3 fixed mesh pockets and repositionable elastic loop panel for secure equipment storage. Available in 500D CORDURA® nylon.

The Rescue Task Force Rapid Response Kit is equipped to provide initial triage of an MCI and provide treatment for 4 or more patients. MOLLE'd on the exterior are 4 C-A-T® tourniquets with holders and nylon dispenser equipped with 4 rolls of color-coded triage tape. Available in Red or Black nylon colors.

### **Features:**

- First Responder Go Bag with modular, personalized organization of mission essential gear AND medical equipment
- Rugged nylon shoulder bag features adjustable, padded strap with quick release buckle
- Multiple storage capabilities including loop-tape/MOLLE attachment system, 2 extra pockets & elastic storage loops with color-coded labels to accommodate up to 8 chem lights (not included)
- Interior compartment features 3 fixed mesh pockets and repositionable elastic loop panel for secure equipment storage
- Also available as a bag only with shoulder bag, 4 TQ holders, and Triage Tape Dispenser (no contents included)

#### **Rapid Response Kit - RTF Edition**

| Color | Item #  | NSN# |
|-------|---------|------|
| BLK   | 85-0716 |      |
| RED   | 85-0717 |      |

#### **Rapid Response Bag only - RTF Edition**

| Color | Item #  | NSN# |
|-------|---------|------|
| BLK   | 80-0412 |      |
| RED   | 80-0413 |      |



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# RAPID RESPONSE KIT *Rescue Task Force edition*

REV100421

## 85-0717 (Red) Kit Contents:



## 85-0716 (Black) Kit Contents:



## Kit Contents:

- 1 x Rapid Response Shoulder Bag
- 4 x Tourniquet Holders
- 1 x Triage Tape Dispenser
- 4 x Color-coded rolls of triage tape
- 4 x C-A-T® Tourniquet
- 4 x Emergency Trauma Dressing, 6 in.
- 4 x NAR Compressed Gauze
- 4 x HyFin® Vent Chest Seal, Twin Pack
- 4 x pair Bear Claw™ Nitrile Gloves, LG
- 1 x NAR Trauma Shears, 7.25 in.

## Dimensions:

- Kit Closed: H 12 in. x W 12 in. x D 6 in.
- Adjustable shoulder strap: 28 in. to 50 in.
- Weight (with contents): 6 lb 4.5 oz



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A   
**MISSION®**