

# NAR BP/STETHOSCOPE COMBO KIT

REV051122

This ultra-compact, multi-pocketed kit contains a set of 4 blood pressure cuffs and stethoscope combo which are instantly accessible and can be stored in most trauma kits. This combo kit features a Stethoscope from ADC with extra large bell and ultra-sensitive diaphragm coupled with a set of four rugged nylon cuffs. The cuff set features adult, large adult, child and infant size cuffs and one-tube bladder with quick release luer slip connector, as well as an ambidextrous palm gauge with high impact casing and trigger valve.

## Specifications:

- SEL#: 09ME-03-STET
- Ultra-compact blood-pressure cuffs in a set of 4 (adult, large adult, child, infant)
- One-tube bladder with quick-release luer slip-connector
- Ambidextrous palm gauge with high-impact casing and trigger valve
- ADC stethoscope with extra-large bell and ultra-sensitive diaphragm

## Kit Contents:

- 4 x Ultra-compact blood pressure cuffs (adult, large adult, child, infant)
- 1 x Ambidextrous palm gauge with high-impact casing and trigger valve
- 1 x One-tube bladder with quick-release luer slip-connector
- 1 x ADC stethoscope with extra-large bell and ultra-sensitive diaphragm
- 1 x Multi-pocketed, heavy-duty, zippered nylon case

## Dimensions:

- Case Closed - L 6.75 in. x W 9 in. x H 5 in.
- Weight: 2 lb



## NAR BP/Stethoscope Combo Kit

ITEM #	NSN #
20-0057	6515-01-542-7723



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**