

URINARY CATHETER KIT

REV041922

The Urinary Catheter Kit contains all components required to utilize a Foley catheter during Prolonged Casualty or Field Care. Each kit comes with Silicone Coated 2-Way Foley catheter, Chloraprep swabsticks, sterile water syringe, 1000ml disposable urinary leg bag, surgical latex gloves, moisture proof underpad, fenestrated drape, and surgical lube.

Specifications:

- Compact kit allowing for ease of use in the field
- Bladder decompression for patients with a bladder obstruction

Dimensions:

- L 12 in. x W 7 in. x D 1.5 in.
- Weight: 7 oz



Suprapubic Catheterization Kit Rx

ITEM #	NSN #
83-0020	--

Kit Contents:



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH  **MISSION®**