

# TOOTH SLOOTH FACTURED TOOTH DETECTOR Pack of 4

REV022522

## Component of the Dental Emergency Response Kit

The Tooth Slooth Fractured Tooth Detector is a component of the Dental Response Aid Kit. This fractured tooth detector is used to pinpoint suspected cusps that cause pain during chewing and a diagnostic device for detection of a cracked tooth or teeth.

### Specifications:

- Diagnostic device for detection of a cracked tooth or teeth
- Pinpoint suspected cusps that cause pain
- Single use and disposable
- Pack of 4

### Dimensions:

- L 6.25 in. x W 6.25 in.
- Weight: 0.07 lbs

### Tooth Slooth Fractured Tooth Detector ☒

ITEM #	QUANTITY	NSN #
ZZ-1415	4 Pack	--



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**