

27 GAUGE DENTAL NEEDLE KIT *Sterile*

REV022522

Component of the Dental Emergency Response Kit

The 27 gauge Dental Needle Kit (5 needles per package) is a component of the Dental Emergency Response Kit. The needle is sterile, disposable, self-threading, and self-aligning universal hub to ensure a consistent fit. Ultrasharp, tribeveled, anticoring lancet point.

Specifications:

- 27 Gauge needle
- Sterile, single use and disposable
- Needle Length: 1 3/8"
- 5 needles per package

Dimensions Each:

- L 9 in. x W 4 in.
- Weight: 0.64 oz

Dimensions 10 Pack:

- L 9 in. x W 4 in.
- Weight: 6.4 oz



27 Gauge Dental Needle Kit **STERILE**

ITEM #	QUANTITY	NSN #
ZZ-1413	5 Needles	--
ZZ-1432	50 Needles	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**