

# DENTAL OPERATIVE KIT *Sterile*

REV070722

## **Component of the Dental Emergency Response Kit**

The Dental Operative Kit contains 3 disposable dental instruments to provide dental care in any location. These tools provide Independent medical or dental healthcare professionals the capabilities to care for dental emergencies at sea, isolated assignments, and austere environments. The Dental Operative Kit contains 1 Rigid Spatula #324 that is used mixing restorative materials and cements, 1 Woodson # 2 Plugger that is used to transport and place dental cements, resins, temporaries and insulating and pulp-capping materials and the Interproximal Carver – IPC carving proximal (in between) tooth surfaces. Used to shape, form, or cut tooth anatomy into amalgam restorations.

### **Specifications:**

- This product requires Medical Device Authorization for purchase
- Sterile, single use and disposable

### **Kit Contents:**

- 1 x Rigid Spatula # 324
- 1 x Woodson # 2 Plugger
- 1 x Interproximal Carver (IPC)

### **Dimensions Each:**

- L 9 in. x W 4 in.
- Weight: 0.14 lbs

### **Dimensions 10 Pack:**

- Weight: 1.4 lbs



### **Dental Operative Kit**

Rx ⓘ **STERILE**

| ITEM #  | QUANTITY | NSN #            |
|---------|----------|------------------|
| ZZ-1411 | Each     | 6520-01-701-5743 |
| ZZ-1430 | 10 Pack  | --               |



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**