

MERET EFAK™ PRO X KITS

REV081021



Compact and hands free kit with a lot of capacity and customizable configuration...

The EFAK™ PRO X EMT first aid kit is designed to help you easily carry and deploy EMT supplies in a very compact, hands free kit, with a lot of capacity and customizable configuration. This kit delivers the tools needed in critical point-of-wounding care for penetrating and other traumatic injuries.

Response ready, this pack provides easy access to gloves and other supplies without opening the main compartment. Gloves can be easily removed from the outer access hole while remaining organized in an inner elastic pouch. Shears, pens, pen lights, radios, pagers, light sticks, cleaning agents, protection sprays, and a rescue knife can be inserted and secured into elastic sleeves on the front, top, and sides of the EFAK™ PRO X. Flip open the main compartment to secure N95 masks, shoe/ head covers, eye protection, hand sanitizer, bandages, and more.

The EFAK™ PRO X™ is available in high viz, blue, red, and black infection control, constructed from waterproof, durable, wipe-clean tarpaulin, and are designed specifically to easily clean and remove blood borne pathogens.

Features:

- Double stitched seams on all stress points
- Exterior hook and loop section for name tape or patches
- Pocket to secure shears
- Internal glove pouch with exterior access
- External elastic sleeves secures additional supplies

Specifications:

- H 4.50 in. x W 6.25 in. x D 1.00 in.
- Weight: 0.95 lbs

Nylon Colors:

BLK

BLU

RED

YEL

Meret EFAK™ Pro X Kits

Color	Level	Item#
BLK	Basic	80-1200
BLU	Basic	80-1201
RED	Basic	80-1202
YEL	Basic	80-1203
BLK	w/ CGLE	80-1204
BLU	w/ CGLE	80-1205
RED	w/ CGLE	80-1206
YEL	w/ CGLE	80-1207



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

Basic



Kit Contents:

- 1 x C-A-T® Tourniquet - Black
- 1 x NAR Wound Packing Gauze
- 1 x 4 in. ETD™ Flat Responder
- 1 x HyFin® Vent Compact, Twin Pack
- 1 x MicroShield CPR Face Mask
- 1 x 6.25 in. Trauma Shears
- 1 x Responder Nitrile Gloves, Lg

Basic with Combat Gauze LE



Kit Contents:

- 1 x C-A-T® Tourniquet - Black
- 1 x QuickClot® Combat Gauze LE
- 1 x 4 in. ETD™ Flat Responder
- 1 x HyFin® Vent Compact, Twin Pack
- 1 x MicroShield CPR Face Mask
- 1 x 6.25 in. Trauma Shears
- 1 x Responder Nitrile Gloves, Lg



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A
MISSION®