

# RESPIRATOR MASK *(Box of 35)*

REV08021



## **A Must Have...**

The FluidShield N95 respirator mask is a fluid-resistant mask approved by the NIOSH that delivers a high level of protection for users, keeping out 95% (0.3 microns) of airborne particles. Exceeding CDC and OSHA specified criteria, this protective mask is durable and lightweight providing a secure seal and 360 degrees of protection. The shapable nose bridge and elastic ear loops provide a custom fit that delivers all-day comfort allowing for longer wear times.

## **Features:**

- NIOSH-approved fluid-resistant mask delivers highest level of protection available
- Required for treating patients with conditions requiring airborne transmission precautions, such as tuberculosis or influenza
- Ideal for exams with known high risk patients
- Procedures involving high levels of spray, splatter and airborne particles

## **Specifications:**

- Dimensions *(box of 35)*: 10 in. x 4.5 in. x 4.25 in.
- Weight *(box of 35)*: 0.75 lb

## **Respirator Mask**

Item #	NSN #
ZZ-1462	



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**