

# NAVAL GUN RESPONSE AID KIT (TRAINER)

REV090122



**Nylon Color:** ODG

## *Slim-line lightweight kit stores easily*

The Naval Gun Response Aid Kit – Trainer is designed to meet the TCCC requirements of the U.S. Navy for the emergency care equipment onboard U.S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels. Ergo dynamic, low profile, casualty response system that features 2 identical removable equipment pouches that attach via hook and loop to a lightweight, durable outer bag, thus meeting the requirements for the U.S. Navy's Naval Gun Response Aid Kit – Trainer. The Kit is slim-line, lightweight kit stores easily on a wall, can be carried over the shoulder or bandolier-style.

### **Features:**

- Constructed of Berry Amendment 500D Cordura solution dyed nylon with Ghillie Tex hardware
- Sand and water resistant
- Self repairing zippers with silent 550 pull cords
- Both outer bag and 2 inner pouches are securable using anti-pilferage seals

### **Specifications:**

- Kit Closed:  
H 16 in. x W 11 in. x D 5.5 in.
- Inner Pouch:  
H 15 in. x W 6 in. x D 4 in.
- Weight: 5 lbs

### **Naval Gun Response Aid Kit - TRAINER**

Color	Item#	NSN#
ODG	85-3205	




**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION**

# NAVAL GUN RESPONSE AID KIT (TRAINER)

REV090122

1 x Naval Gun  
Response Bag



2 x Internal Pouches  
(Naval First Aid Box  
Response Kit)

## Kit Contents:

- 2 x Naval First Aid Box Response Pouch
- 2 x C-A-T® Tourniquet, Trainer
- 2 x HyFin® Vent Chest Seal Trainer, Twin Pack
- 4 x ETD™ 6 in. Emergency Trauma Dressing
- 2 x NAR Compressed Gauze
- 4 x Elastic Wrap Bandage (3 in. x 5 yd.)
- 2 x NAR Wound Packing Gauze
- 2 x Triangular Bandage Kit - ODG
- 2 x ARS® Needle, 14 ga.
- 2 x Nasopharyngeal Airway 28F
- 2 x BurnTec Dressing (4 in. x 16 in.)
- 2 x SAM® Splint II
- 2 x Microshield CPR Mask
- 2 x Fox Eye Shield with Garter
- 2 x pair Bear Claw Nitrile Trauma Gloves
- 2 x Trauma Shears, Large
- 2 x Surgical Lubricant (Single Use)
- 2 x 2 in. Surgical Tape
- 2 x TCCC Triage Card
- 2 x ChemLight™, Green
- 2 x NAR Survival Blanket
- 2 x Permanent Marker, Large



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**