

# OSL™ OVERSIZED LITTER *Specialty Litter*

REV070822

## ***Specialty litter designed for safely evacuating larger patients...***

Designed to safely evacuate larger patients with both comfort and ease. This design has a 30% increase in patient surface area compared to standard-sized evacuation platforms, extra padded carrying handles, flame-retardant cover, and additional stirrups offering firmer patient support as well as better weight distribution on the floor.

Like the Raven 90C, the OSL has a bi-fold configuration. The cover material is a very durable Nomex/Kevlar blend.

Also, a supplemental Back Rest (not shown) to support the patient at both 90° and 45° angles may be purchased separately.

### **Specifications:**

- Folded: H 45 in. x W 9 in. x D 8 in.
- Open: L 90 in. x W 31 in. x D 6 in.
- Weight: 22 lb 6 oz (10.2 kg)
- Static Payload: 2000 lb (909 kg)
- Ground Clearance: 1.5 in.

### **Features:**

- Certified safe-to-fly by the U.S. Air Force
- Flame-retardant cover made of a durable Nomex/Kevlar-blend fabric
- 6 stirrups for firmer patient support and better weight distribution
- Ergonomically designed handles
- Additional carrying loops
- 6 IV attachment points
- 3 patient securing straps
- Attached storage pouch, Bifold configuration for easy storage
- U.S. Patent 5,598,592, Patent(s) Pending



*Designed with 6 stirrups for firmer patient support and better weight distribution*

### ***OSL OverSized Litter***

Item#	NSN#
60-0023	6530-01-552-7708



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**