

# MEDIC TRAUMA PACK BAG *Bag Only*

REV103123



## **Tactical Healthcare Professional Medic Trauma Bag**

The Medic/Trauma Pack BAG ONLY is a light framed platform with silent, dirt and sand-resistant zippers that close over a large main compartment with modular, removable pull bags for all your trauma management essentials.

The adjustable shoulder and sternum straps make for a comfortable fit. Cinch straps are included for the sides to maintain load integrity during stress maneuvers.

Grommets are provided at the base of the pack for drainage.

### **Features:**

- Over 1000 cubic inches of gear capacity
- Contoured, padded shoulder straps for added comfort
- Constructed of durable 500 Denier Nylon
- Sand & water resistant zipper with silent cord pulls
- External Cinch Straps to tighten pack and decrease load shift
- MOLLE Style Attachment System Compatible

### **Bag Specs:**

- H 19 in. x W 19 in. x D 10 in.

#### **Medic Trauma Pack Bag**

ITEM #	NSN#
80-0073W	
80-0075W	----

#### **Available Colors:**



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**