

BRAIDED ROLL DENTAL SPONGE *Non-Sterile*

REV031022

Component of the Dental Emergency Response Kit

The Dental Rolled Braided Sponge is non-sterile .375" X 1.5" and is a component of the Dental Emergency Response Kit. A braided wrapped cotton roll medium non-sterile and disposable.

Specifications:

- 50 non-sterile individual rolls per pack
- Starch-free, so they will not stick to the mucosa
- 100% cotton interior wrapped in an absorbent non-woven fabric

Dimensions:

- L 1.5 in. x W 0.375 in.
- Weight: 0.11 lbs



Braided Roll Dental Sponge

ITEM #	QUANTITY	NSN #
ZZ-1414	50 Rolls	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**