

CRO MEDICAL PELVIC BINDER

REV031022

Reduce and Stabilize Open-Book Pelvic Fractures

Easily reduce and stabilize open-book pelvic fractures using the patented tightening system. The CRO pelvic binder packs flat and is designed to fit 95% of the population using the Standard size.

BOA precision dial system has proven reliability and testing to easily achieve at least 33lbs (150N) of circumferential force.

The BOA hardware was selected to eliminate human error during application. Here's how: The average human hand can tighten between 40-50 lbs. There are no cases of over reduction in the CPG's, so achieving pressure over the minimum is most desired. How do you know you have enough pressure? There's no thought involved, just tighten until you can't tighten anymore.

The CoTCCC update, released on 31JAN2017 indicates pelvic binding under "Circulation" for bleeding control, due to common high-energy battlefield MOI's associated with significant morbidity and mortality due to associated vascular injuries. 26% of service members that died in OEF/OIF died with a pelvic fracture.

Specifications:

- BOA Precision Dial System
- Hook and loop retention for Pressure Delivery Device
- Approximately half the size in cubic inches of other binders
- Lightweight 4 in. double-back wrap
- Easy, one-hand tightening
- Low-light, no-light features

Dimensions:

- Standard Waist Size: 32 in. - 50 in. hips
- Weight: 0.75 lbs



CRO Medical Pelvic Binder

ITEM #	QUANTITY	NSN #
50-0044	1 Each	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®