

DENTAL SCALER KIT *Sterile*

REV092722

Component of the Dental Emergency Response Kit

The Dental Scaler Kit contains two disposable dental instruments to provide Dental Care for independent medical or dental healthcare professionals the capabilities to care for dental emergencies at sea, isolated assignments, and austere environments. The Dental Scaler Kit contains 1 15/30 Towner designed for the removal of supragingival calculus and 1 McCall Curette 13/14 designed with elliptical blade shape with a rounded toe to remove both supragingival and subgingival calculus.

Specifications:

- Sterile, single use and disposable
- Ideal for dental emergencies in austere environments

Kit Contents:


- 1 X 15/30 Towner
- 1 X McCall Curette

Dimensions Each:

- Packaged - L 7.75 in. x W 3.75 in.
- Weight: 0.12 lbs

Dimensions 10 Pack:

- Packaged - L 19.5 in. x W 15.5 in. x H 14 in.
- Weight: 1.1 lbs



Dental Scaler Kit

Rx  **STERILE**

ITEM #	QUANTITY	NSN #
ZZ-1410	Each	6520-01-701-6776
ZZ-1429	10 Pack	--



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**