

QUIKCLOT EMS 4X4 HEMOSTATIC GAUZE

REV012023



Effective, simple and safe when trying to control life-threatening hemorrhage...

QuikClot® EMS 4 x 4 dressing is a soft, white, sterile, non woven 4 in. by 4 in. 4-ply gauze impregnated with the hemostatic, kaolin, an inert mineral that does not generate heat.

Each QuikClot EMS 4 X 4 dressing is individually wrapped in an easy rip, pouch. Indicated for temporary external control of traumatic bleeding, QuikClot EMS 4 X 4 dressing is flexible and pliable and contours to all wounds.

Effective, simple and safe to use when trying to control life-threatening hemorrhage from wounds that are not amenable to tourniquet application.

Features:

- **Controls Bleeding:** Promotes clotting within minutes, allowing responders to control bleeding in an efficient and timely manner.
- **Easy To Use:** QuikClot EMS 4 x 4 Hemostatic Gauze uses the same placement technique as standard gauze.
- **Heat Free:** QuikClot products produce NO Heat
- **Dimensions:**
 - Packaged: H 5.5 in. x W 5 in. x D 0.2 in.
 - Gauze Dimensions: W 4 in.x L 4 in.
 - Weight: 0.8 oz

QuikClot Combat Gauze LE

Item #	NSN#
30-0095	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A
MISSION®