

Medical Protective **FACE SHIELDS**

REV090120

Designed to work in conjunction with a face mask to provide another level of PPE...

The Medical Protective Face Shield is designed as a cost effective general splash guard to add an extra layer of personal protection from debris and chemical hazards. This professional face shield is simply designed to work in conjunction with a face mask to provide another level of PPE. The OSHA compliant clear shield is constructed of 14 mil BPA free plastic with a fully adjustable, one size fits all headband that can accommodate extra room between the face and the shield for goggles. Made in the USA.

Features:

- Cost effective, simply designed professional face shield to protect the face and eyes from flying debris and chemical hazards
- Adjustable headband allows leaving room for goggles by holding the shield away from the face
- Available in convenient box of 20
- Ideal for pandemic preparedness stockpiling
- Made in the USA

Specifications:

- Shield: H 8.5 in. x W 13.5 in.
- Black Strap Headband: W 1.25 in. x L 26 in.
- Weight: 1.36 oz

Protective Face Shields

ITEM#	UNIT
70-0262	Box of 20

Face Mask & Face Shield are sold separately



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH

MISSION®