

# KN95 RESPIRATOR MASK (Box of 50)

REV091520



## Specifications:

- Dimensions: 10 in. x 4.75 in. x 4.75 in. (box of 50)
- Weight: 0.32 oz (each mask)

## A Must Have...

Offering 5-filtration-layers which meet and exceed KN95 standards for protection, our KN95 Respirator face mask will keep out 95% (0.3 microns) of airborne particles for maximum protection and perform better than a standard disposable 3-Ply mask or conventional face cover. This protective face mask is lightweight and comfortable, with ear loops and a shapeable nose bridge for a flexible and secure seal. Features include 360° protection with three dimensional ventilation space for more comfortable breathing.

## Features:

- 5-Layers of protection blocking at least 95% (0.3 microns) of airborne particulates
- More protection than standard disposable 3-Ply masks or your conventional face cover
- Lightweight and foldable, easy to store and carry
- 360° protection with three dimensional ventilation space for more comfortable breathing
- Equipped with metal fastening strip ensuring optimal sealing between the bridge of the nose and the mask
- Sold as a box of 50

## KN95 Respirator Mask

Item #	NSN #
ZZ-1308	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**