

PROTECTIVE ISOLATION GOWN

REV082120



A Must Have...

North American Rescue is proud to offer the Protective Isolation Gown. It meets FDA Level 1, 2, and 3 protective gown guidance to help mitigate the spread of communicable diseases. They are individually wrapped for cleanliness and feature front or rear break-away capability to remove.

Features:

- Evaluated and passed:
 - ASTM F2407-06
 - Liquid Penetration Test ATTCC 42 (*American Association of Textile Chemists and Colorists*)
- Meets FDA Level 1, 2, and 3 protective gown guidance
- Front or rear break-away capability to remove
- Individually wrapped for cleanliness
- 1.25 mil LDPE

Specifications:

- Dimensions: L 49 in. x W 36 in. (*at the shoulders*)
- Dimensions (*Arms*): L 20 in. x W 10.5 in.

Protective Isolation Gown

Item #	NSN #
70-0263	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®