

IMPERVIOUS ISOLATION GOWNS *Pack of 10*

REV081720



Protection against splashing and spraying of fluids...

Made with single-layer chlorinated polyethylene (CPE), a thermoplastic material. Lightweight, durable, and resistant to flames and chemicals. Knee-length blue gowns feature long-sleeves, cuff bands with thumb hooks, and waist ties. Pullover design. One size fits most. Not made with natural rubber latex. Sold in a pack of 10.

Features:

- Knee-length pullover gown features long-sleeves, cuff bands with thumb hooks, and waist ties
- Provides protection against splashing and spraying of body fluids
- Impervious to fluids with single-layer thermoplastic material - chlorinated polyethylene (CPE)
- Lightweight, durable, and resistant to flames and chemicals
- One size fits most, Non Sterile, Latex Free, Single Use, Blue

Specifications:

- Gown (at shoulders) : L 45 in. x W 29 in.
- Arms: L 23 in. x W 12 in.
- Packaged (10 pack) : L 12 in. x W 11 in.
- Weight (10 pack) : 1.17 lb

Impervious Isolation Gowns (Pack of 10)

Color	Item#	NSN#
Blue	ZZ-0750	



***Featuring convenient
Thumb Hook Loops***



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**